Instructions for use

- EN** see www.fkg.ch for additional information
FR informations supplémentaires sur www.fkg.ch
DE weitere Informationen auf www.fkg.ch

Golden rules


- Speed: 600 rpm
 Torque: 1.5 Ncm
- Long back and forth gentle strokes
 - Light touch, let the instrument work
 - Work 3-4 seconds in a row, get out, clean the blade and irrigate the canal

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iRace 1 blister includes/contient/enthält: 1 × R1, 1 × R2, 1 × R3

For most cases

- straight, slightly curved and/or large
Prepare glide path first, then

1. Introduce R1 in rotation and reach working length (WL); In case R1 does not reach the WL, do not force and go to step 1 of iRace Plus protocol
2. Continue shaping with R2 up to WL;
3. Finish the shaping with R3 up to WL.

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iRace Plus 1 blister includes/contient/enthält: 2 × (R1a + R1b)

For difficult cases

- severely curved, narrow and/or calcified
Following of iRace protocol step 1

1. Use R1a to reach WL;
2. Continue shaping with R1b up to WL;

When WL is reached, continue from step 2 of iRace protocol.

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iRace Plus

SafetyMemoDisc (SMD)

Pull off petals to master metal fatigue:

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- 21/16mm ● 25/16mm ● 31/16mm

EN Sterilization protocol for reusable medical devices

Non sterile products

First use

Carefully follow the rules herebelow

Cleaning and Sterilization

Cleaning

1. Manual cleaning with/without Ultrasonic assistance

Thorough cleaning allows efficient disinfection & sterilization.

2. Rinse

Under distilled/demineralized water at least 1 min.

Then dry the instruments.

Sterilization

5. According EN/ISO 17664 protocol

Autoclave: 134°C/273°F,

2.2 bar during 18 min.

 Always refer to manufacturer's instructions for cycle & duration

Second and following uses

7. Pre-disinfection

Plunge instruments after use in detergent solution and eventually brush them manually

 Complete operating protocol and warnings on www.fkg.ch

General information

Pre-/Disinfection

DO NOT USE

solutions containing:

- phenol (corrosion)
- aldehyde (blood fixation)
- di-/trietanolamines (corrosion)

Caution:

This product contains nickel and should not be used for individuals with known allergic sensitivity to this metal.

A single patient use is recommended to avoid cross-contamination.

Sterile products

First use

1. Open the alveolus
2. Extract the instrument (use of gloves imperative)
3. Use the instrument according to good dental practice

Second use:

Carefully follow the rules herebelow

Second and following uses.

3. Inspection - Check

Sort out damaged or worn out instruments.

4. Packing

Place the instruments in a suitable support/container and pack the devices in sterilisation pouches ISO 11607-1.

6. Storage

Keep devices in sterilization wrap/pouch in a dry and clean environment.

8. Rinse

Under running water at least 1 min.

Then dry the instruments.

9. Follow steps 1 - 6 above

Sterilization with chemiclav or dry air devices

have not been validated by FKG Dentaire.

Refer to manufacturer's instructions for cycle & duration. Nevertheless our instruments withstand such methods.

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