



**EXPAND  
YOUR  
MIND**

[www.fkg.ch/xpendo](http://www.fkg.ch/xpendo)

**SCIENTIFIC OVERVIEW**  
12/2019

MECHANICAL PROPERTIES	DEBRIS/SMEAR LAYER REMOVAL	CLINICAL TRIAL/CASE	BACTERIA REMOVAL	RETREATMENT	APICAL DEBRIS EXTRUSION	TESTED AT BODY TEMPERATURE	
<b>2019</b>							
			●			✓	<p style="text-align: center;"><b>(1) Teves A, et al. JCED. 2019</b></p> <p><b>Compared with:</b> positive pressure irrigation (PPI).  <b>Comments:</b> Biofilm (<i>Enterococcus faecalis</i>, <i>Eikenella Corrodens</i> and <i>Streptococcus anginosus</i>) removal was more efficient with NaOCl 4% than with chlorhexidine 2%. In addition, the use of XP-endo Finisher improved the biofilm removal efficiency further.</p>
	●					✓	<p style="text-align: center;"><b>(2) Marques-da-Silva B, et al. IEJ. 2019</b></p> <p><b>Compared with:</b> passive ultrasonic irrigation (PUI), EDDY and EndoActivator.  <b>Comments:</b> XP-endo Shaper &amp; Finisher and EDDY showed better results than the other groups in regards on calcium hydroxide removal.</p>
	●					✓	<p style="text-align: center;"><b>(3) De-Deus G, et al. COI. 2019</b></p> <p><b>Compared with:</b> passive ultrasonic irrigation (PUI).  <b>Comments:</b> Overall, XP-endo Finisher and PUI were equally efficient in debris removal.</p>
				●		✓	<p style="text-align: center;"><b>(4) Aksel H, et al. IEJ. 2019</b></p> <p><b>Compared with:</b> none.  <b>Comments:</b> The use of XP-endo Finisher significantly reduced the remaining volume of filling material after the initial retreatment procedure with Protaper Universal Retreatment.</p>
			●			✓	<p style="text-align: center;"><b>(5) Pacheco-Yanes J, et al. COI. 2019</b></p> <p><b>Compared with:</b> non-agitated irrigation and passive ultrasonic irrigation (PUI).  <b>Comments:</b> The XP-endo Finisher instrument highlighted a significantly better distribution of the irrigant when compared to PUI and control.</p>

● All comparisons are at least equal to XP-endo Finisher  
 ● At least 1 comparison is negative for XP-endo Finisher  
 ● All comparisons are negative for XP-endo Finisher

✓ Yes    ✗ No

MECHANICAL	DEBRIS/SMEAR LAYER	CLINICAL	BACTERIA	RETREATMENT	EXTRUSION	BODY TEMPERATURE	
			●				<p>(6) <b>De-Deus G, et al. IEJ. 2019</b></p> <p><b>Compared with:</b> passive ultrasonic irrigation (PUI).  <b>Comments:</b> Significantly more root filling material was removed with XP-endo Finisher R (32%) than with PUI (12%).</p>
	●						<p>(7) <b>Azimian S, et al. DRJ. 2019</b></p> <p><b>Compared with:</b> none.  <b>Comments:</b> The authors claimed that XP-endo Finisher has no superiority compared to control group. However, XP-endo Finisher was not used in combination with irrigant nor were the tests performed at body temperature.</p>
	●						<p>(8) <b>Jayakumaar A, et al. IJDR. 2019</b></p> <p><b>Compared with:</b> conventional irrigation.  <b>Comments:</b> The debris and smear layer scores were significantly in favor of using XP-endo Finisher after instrumentation with Hyflex and Protaper Next.</p>
			●				<p>(9) <b>Campello AF, et al. IEJ. 2019</b></p> <p><b>Compared with:</b> solvent.  <b>Comments:</b> Solvent was not of additional help, but XP-endo Finisher R removed a significant amount of (additional) gutta percha after Mtwo instrumentation.</p>
			●				<p>(10) <b>Machado AG, et al. IEJ. 2019</b></p> <p><b>Compared with:</b> none.  <b>Comments:</b> As for the number of cases with total filling material removal, XP-endo Shaper was associated with better results than TRUShape. The supplementary approach with the XP-endo Finisher R instrument significantly improved removal of existing filling material.</p>

MECHANICAL	DEBRIS/SMEAR LAYER	CLINICAL	BACTERIA	RETREATMENT	EXTRUSION	BODY TEMPERATURE	
			●				<p>(11) Zhao Y, et al. IEJ. 2019</p> <p><b>Compared with:</b> conventional irrigation and passive ultrasonic irrigation (PUI).</p> <p><b>Comments:</b> The additional irrigation (after instrumentation with XP-endo Shaper or Reciproc Blue) showed significantly less remaining dentinal debris. Conventional irrigation showed less reduction than PUI and XP-endo Finisher when used after Reciproc Blue. PUI and Xp-endo Finisher showed no difference.</p>
			●				<p>(12) Sasanakul P, et al. JOE. 2019</p> <p><b>Compared with:</b> standard manual instrumentation (CF), minimal manual instrumentation (MI), Navitip FX needle (NFX), Non-agitated irrigation (NI), Passive ultrasonic irrigation (PUI) and SAF.</p> <p><b>Comments:</b> The bacteria removal hierarchy was as follow: MI &gt; NFX = XP-endo Finisher &gt; CF = SAF = PUI &gt; NI. The samples were not harmonized in instrumentation sizes/conditions prior the irrigation/disinfection. There is no indication in the article that the experiment was performed at body temperature.</p>
			●				<p>(13) Carvalho MC, et al. BOR. 2019</p> <p><b>Compared with:</b> none.</p> <p><b>Comments:</b> XP-endo Finisher significantly reduced the bacterial load after instrumentation with XP-endo Shaper or Reciproc Blue.</p>

2018						
MECHANICAL	DEBRIS/SMEAR LAYER	CLINICAL	BACTERIA	RETREATMENT	EXTRUSION	BODY TEMPERATURE
		●				<p>(14) <b>Rodrigues E-A, et al. JCED. 2018</b></p> <p><b>Compared with:</b> none; case report.  <b>Comments:</b> A 22-year-old patient was treated for a dens invaginatus type II. The chemo-mechanical treatment was achieved with XP-endo Finisher and NaOCl. At the 14 months follow-up, the patient was asymptomatic and osseous healing of the lesion was visible on the X-ray.</p>
	●					<p>(15) <b>Keskin C, et al. JDRDCDP. 2018</b></p> <p><b>Compared with:</b> conventional irrigation and passive ultrasonic irrigation (PUI).  <b>Comments:</b> XP-endo Finisher and PUI showed significantly better triple antibiotic paste (TAP) removal than conventional needle irrigation at all time points. At 7- and 21-days TAP incubation time, XP-endo Finisher and PUI presented similar results. Only at day 90, PUI removed significantly more TAP than XP-endo Finisher. There is no indication that the study was performed at body temperature.</p>
	●					<p>(16) <b>Ulusoy Öi, et al. IEJ. 2018</b></p> <p><b>Compared with:</b> non-agitated irrigation and passive ultrasonic irrigation (PUI).  <b>Comments:</b> XP-endo Finisher was the most effective activation method to remove organic tissue. The irrigation solutions (NaOCl, NaOCl-EDTA, and NaOCl + HEBP) did not significantly influence the results.</p>

MECHANICAL	DEBRIS/SMEAR LAYER	CLINICAL	BACTERIA	RETREATMENT	EXTRUSION	BODY TEMPERATURE	
			●				<p>(17) <b>Bedier MM, et al. RDE. 2018</b></p> <p><b>Compared with:</b> conventional irrigation.  <b>Comments:</b> Effectiveness of iRace or XP-endo Shaper + XP-endo Finisher to remove artificially grown <i>E. faecalis</i> from the middle part of the canal and tubules compared to the same initial files and conventional irrigation only. All groups showed superior results for XP-endo Finisher compared to conventional irrigation.</p>
	●						<p>(18) <b>Kfir A, et al. AEJ. 2018</b></p> <p><b>Compared with:</b> conventional irrigation, passive ultrasonic irrigation (PUI) and SAF.  <b>Comments:</b> XP-endo Finisher, PUI, and SAF removed similar amount of <math>\text{Ca(OH)}_2</math> from artificial grooves. Conventional removed significantly less <math>\text{Ca(OH)}_2</math> than the other groups.</p>
●							<p>(19) <b>Vaz-Garcia ES, et al. BDJ. 2018</b></p> <p><b>Compared with:</b> XP-Clean.  <b>Comments:</b> XP-endo Finisher instruments showed improved performance when compared with XP-Clean instruments, demonstrating higher cyclic fatigue resistance and lower roughness.</p>
					●		<p>(20) <b>Azim AA, et al. COI. 2018</b></p> <p><b>Compared with:</b> conventional irrigation, EndoActivator, EndoVac and PIPS.  <b>Comments:</b> Only negative pressure (EndoVac) had no extrusion. The amount of extruded debris, for the other methods, was similar. This experimental setup has low clinical significance since the pressure of the periapical tissues are not taken into account.</p>



							MECHANICAL	DEBRIS/SMEAR LAYER	CLINICAL	BACTERIA	RETREATMENT	EXTRUSION	BODY TEMPERATURE
													<p>(21) Kfir A, et al. COI. 2018</p> <p><b>Compared with:</b> SAF.</p> <p><b>Comments:</b> The complete SAF sequence extruded less than the combination of Proglider, Protaper Next, and XP-endo Finisher. No significant difference was observed between Proglider and XP-endo Finisher. No statistics between SAF and XP-endo Finisher was provided.</p>
													<p>(22) Silva EJNL, et al. IEJ. 2018</p> <p><b>Compared with:</b> none.</p> <p><b>Comments:</b> XP-endo Finisher and XP-endo Finisher R significantly removed filling material after initial retreatment with round instruments. No significant difference was observed between XP-endo Finisher and XP-endo Finisher R.</p>
<b>2017</b>													
													<p>(23) Zand V, et al. JCED. 2017</p> <p><b>Compared with:</b> none.</p> <p><b>Comments:</b> For smear layer removal, XP-endo Finisher with 17% EDTA for one minute was the most effective combination.</p>
													<p>(24) Turkydin D, et al. JOE. 2017</p> <p><b>Compared with:</b> conventional irrigation and passive ultrasonic irrigation (PUI).</p> <p><b>Comments:</b> XP-endo Finisher showed significantly lower amount of remaining triple antibiotic past (TAP) than needle irrigation and PUI.</p>

	MECHANICAL	DEBRIS/SMEAR LAYER	CLINICAL	BACTERIA	RETREATMENT	EXTRUSION	BODY TEMPERATURE
	●					✗	<p>(25) <b>Uygun AD, et al. AEJ. 2017</b></p> <p><b>Compared with:</b> conventional irrigation, passive ultrasonic irrigation (PUI), and TRUShape.</p> <p><b>Comments:</b> Conventional irrigation had the poorest scores in terms of Ca(OH)<sub>2</sub> removal. XP-endo Finisher, TRUShape, and PUI groups had similar results in removing calcium hydroxide.</p>
	●					✗	<p>(26) <b>Wigler R, et al. IEJ. 2017</b></p> <p><b>Compared with:</b> conventional irrigation and passive ultrasonic irrigation (PUI).</p> <p><b>Comments:</b> XP-endo Finisher and PUI removed significantly more Ca(OH)<sub>2</sub> from artificial grooves than conventional irrigation. There is no indication in the article that the experiment was performed at body temperature.</p>
	●					✓	<p>(27) <b>Hamdan R, et al. JCED. 2017</b></p> <p><b>Compared with:</b> passive ultrasonic irrigation (PUI).</p> <p><b>Comments:</b> XP-endo Finisher showed a significant superiority over PUI in removing Ca(OH)<sub>2</sub> from the apical third after 3 minutes of activation.</p>
	●					✗	<p>(28) <b>Gokturk H, et al. JAOSR. 2017</b></p> <p><b>Compared with:</b> CanalBrush, conventional irrigation, double side vented needle, laser-activated irrigation (LAI), passive ultrasonic irrigation (PUI), and Vibringe.</p> <p><b>Comments:</b> LAI and PUI showed the highest mean rate of Ca(OH)<sub>2</sub> removal from artificial grooves. However, there is no indication that the experiments were performed at body temperature.</p>



	MECHANICAL	DEBRIS/SMEAR LAYER	CLINICAL	BACTERIA	RETREATMENT	EXTRUSION	BODY TEMPERATURE
	●					×	<p>(29) Leoni GB, et al. IEJ. 2017</p> <p><b>Compared with:</b> apical positive pressure (APP), passive ultrasonic irrigation (PUI), and SAF.</p> <p><b>Comments:</b> XP-endo finisher and PUI was equally effective overall. However, XP-endo Finisher was best in the apical area that is the most challenging area for debris removal and disinfection.</p>
				●		×	<p>(30) Karamifar K, et al. IranEJ. 2017</p> <p><b>Compared with:</b> manual instrumentation and RACE.</p> <p><b>Comments:</b> The use of XP-endo Finisher file resulted in cleaner canal walls and was more effective in removing gutta-percha from the coronal toward the apical part of the canal. Emphasizes the advantage of the XP-endo Finisher in the apical part of the canal especially.</p>
	●					×	<p>(31) Elnaghy AM, et al. Odont. 2017</p> <p><b>Compared with:</b> conventional irrigation, EndoActivator, and non-agitated irrigation.</p> <p><b>Comments:</b> XP-endo Finisher and EndoActivator were superior to the other methods tested. There is no indication that the experiment was performed at body temperature.</p>
			●			✓	<p>(32) Bao P, et al. JOE. 2017</p> <p><b>Compared with:</b> conventional irrigation (CNI) and passive ultrasonic irrigation (PUI).</p> <p><b>Comments:</b> The best biofilm removal – inside and outside of artificial grooves – was achieved by XP-endo Finisher. PUI and CNI followed.</p>

MECHANICAL	DEBRIS/SMEAR LAYER	CLINICAL	BACTERIA	RETREATMENT	EXTRUSION	BODY TEMPERATURE	
	●						<p>(33) Keskin C, et al. JOE. 2017</p> <p><b>Compared with:</b> CanalBrush, conventional irrigation (SI), EndoActivator (EA), and passive ultrasonic irrigation (PUI).  <b>Comments:</b> XP-endo Finisher and PUI removed similar amount of Calcium Hydroxide. Both removed significantly more Calcium Hydroxide than SI, EA, and Canal.</p>
<b>2016</b>							
	●						<p>(34) Gokturk H, et al. JAOSR. 2016</p> <p><b>Compared with:</b> CanalBrush, conventional irrigation, double side vented needle, laser-activated irrigation (LAI), passive ultrasonic irrigation (PUI), and Vibringe.  <b>Comments:</b> The removal of double antibiotic paste from artificial grooves was investigated for various irrigation protocols. Depending on the location into the canal, the significance levels between the protocols differed. There is no indication that the study was performed at body temperature.</p>
			●				<p>(35) Alves FRF, et al. JOE. 2016</p> <p><b>Compared with:</b> passive ultrasonic irrigation (PUI).  <b>Comments:</b> Both XP-endo Finisher and PUI exhibited antibacterial effectiveness, but only the XP-endo Finisher caused a significant reduction in the bacterial counts after chemomechanical preparation.</p>
				●			<p>(36) Alves FRF, et al. JOE. 2016</p> <p><b>Compared with:</b> none.  <b>Comments:</b> XP-endo Finisher removed significantly more debris after initial retreatment with round files.</p>

MECHANICAL	DEBRIS/SMEAR LAYER	CLINICAL	BACTERIA	RETREATMENT	EXTRUSION	BODY TEMPERATURE
			●		✗	<p>(37) Azim AA, et al. JOE. 2016</p> <p><b>Compared with:</b> conventional irrigation, EndoActivator, and PIPS.</p> <p><b>Comments:</b> XP-endo Finisher showed significantly more bacterial reduction (from canal wall to 50 µm inside dentinal tubule) than all other activation methods. Deeper into the tubules, PIPS killed significantly more bacteria.</p>

## References

- Teves A, Blanco D, Casaretto M, Torres J, Alvarado D, Jaramillo DE. Effectiveness of different disinfection techniques of the root canal in the elimination of a multi-species biofilm. *J Clin Exp Dent*. 2019 Nov 1;11(11):e978–83.
- Marques-da-Silva B, Alberton CS, Tomazinho FSF, Gabardo MCL, Duarte M a. H, Vivian RR, et al. Effectiveness of five instruments when removing calcium hydroxide paste from simulated internal root resorption cavities in extracted maxillary central incisors. *Int Endod J*. [Epub ahead of print].
- De-Deus G, Belladonna FG, de Siqueira Zuolo A, Perez R, Carvalho MS, Souza EM, et al. Micro-CT comparison of XP-endo Finisher and passive ultrasonic irrigation as final irrigation protocols on the removal of accumulated hard-tissue debris from oval shaped-canals. *Clin Oral Investig*. 2019 Jul;23(7):3087–93.
- Aksel H, Küçükaya Eren S, Askerbeyli Örs S, Serper A, Ocak M, Çelik HH. Micro-CT evaluation of the removal of root fillings using the ProTaper Universal Retreatment system supplemented by the XP-Endo Finisher file. *Int Endod J*. 2019 Jul;52(7):1070–6.
- Pacheco-Yanes J, Provenzano JC, Marceliano-Alves MF, Gazzaneo I, Pérez AR, Gonçalves LS, et al. Distribution of sodium hypochlorite throughout the mesial root canal system of mandibular molars after adjunctive irrigant activation procedures: a micro-computed tomographic study. *Clin Oral Investig*. 2019 Jun 26;
- De-Deus G, Belladonna FG, Zuolo AS, Cavalcante DM, Carvalho JCA, Simões-Carvalho M, et al. XP-endo Finisher R instrument optimizes the removal of root filling remnants in oval-shaped canals. *Int Endod J*. 2019 Jun;52(6):899–907.
- Azimian S, Bakhtiar H, Azimi S, Esnaashari E. In vitro effect of XP-Endo finisher on the amount of residual debris and smear layer on the root canal walls. *Dent Res J*. 2019 Jun;16(3):179–84.
- Jayakumaar A, Ganesh A, Kalaiselvam R, Rajan M, Deivanayagam K. Evaluation of debris and smear layer removal with XP-endo finisher: A scanning electron microscopic study. *Indian J Dent Res Off Publ Indian Soc Dent Res*. 2019 Jun;30(3):420–3.
- Campello AF, Almeida BM, Franzoni MA, Alves FRF, Marceliano-Alves MF, Rôças IN, et al. Influence of solvent and a supplementary step with a finishing instrument on filling material removal from canals connected by an isthmus. *Int Endod J*. 2019 May;52(5):716–24.
- Machado AG, Guilherme BPS, Provenzano JC, Marceliano-Alves MF, Gonçalves LS, Siqueira JF, et al. Effects of preparation with the Self-Adjusting File, TRUShape and XP-endo Shaper systems, and a supplementary step with XP-endo Finisher R on filling material removal during retreatment of mandibular molar canals. *Int Endod J*. 2019 May;52(5):709–15.
- Zhao Y, Fan W, Xu T, Tay FR, Gutmann JL, Fan B. Evaluation of several instrumentation techniques and irrigation methods on the percentage of untouched canal wall and accumulated dentine debris in C-shaped canals. *Int Endod J*. 2019 Mar 21;

12. Sasanakul P, Ampornaramveth RS, Chivatxaranukul P. Influence of Adjuncts to Irrigation in the Disinfection of Large Root Canals. *J Endod.* 2019 Mar;45(3):332–7.
13. Carvalho MC, Zuolo ML, Arruda-Vasconcelos R, Marinho ACS, Louzada LM, Francisco PA, et al. Effectiveness of XP-Endo Finisher in the reduction of bacterial load in oval-shaped root canals. *Braz Oral Res.* 2019;33:e021.
14. Rodrigues E-A, Belladonna F-G, De-Deus G, Silva E-J-N-L. Endodontic management of type II dens invaginatus with open apex and large periradicular lesion using the XP-endo Finisher: A case report. *J Clin Exp Dent.* 2018 Oct;10(10):e1040–4.
15. Keskin C, Güler DH, Sarıyılmaz E. Effect of intracanal time of triple antibiotic paste on its removal from simulated immature roots using passive ultrasonic irrigation and XP-endo Finisher. *J Dent Res Dent Clin Dent Prospects.* 2018;12(4):288–93.
16. Ulusoy Ö, Savur İG, Alaçam T, Çelik B. The effectiveness of various irrigation protocols on organic tissue removal from simulated internal resorption defects. *Int Endod J.* 2018;51(9):1030–6.
17. Bedier MM, Hashem AAR, Hassan YM. Improved dentin disinfection by combining different-geometry rotary nickel-titanium files in preparing root canals. *Restor Dent Endod.* 2018 Aug 22;43(4).
18. Kfir A, Blau-Venezia N, Goldberger T, Abramovitz I, Wigler R. Efficacy of self-adjusting file, XP-endo finisher and passive ultrasonic irrigation on the removal of calcium hydroxide paste from an artificial standardized groove. *Aust Endod J J Aust Soc Endodontology Inc.* 2018 Apr;44(1):26–31.
19. Vaz-Garcia ES, Vieira VTL, Petitet NP da SF, Moreira EJJ, Lopes HP, Elias CN, et al. Mechanical Properties of Anatomic Finishing Files: XP-Endo Finisher and XP-Clean. *Braz Dent J.* 2018 Mar;29(2):208–13.
20. Azim AA, Aksel H, Margaret Jefferson M, Huang GT-J. Comparison of sodium hypochlorite extrusion by five irrigation systems using an artificial root socket model and a quantitative chemical method. *Clin Oral Investig.* 2018 Mar;22(2):1055–61.
21. Kfir A, Moza-Levi R, Herteanu M, Weissman A, Wigler R. Apical extrusion of debris during the preparation of oval root canals: a comparative study between a full-sequence SAF system and a rotary file system supplemented by XP-endo finisher file. *Clin Oral Investig.* 2018 Mar;22(2):707–13.
22. Silva EJNL, Belladonna FG, Zuolo AS, Rodrigues E, Ehrhardt IC, Souza EM, et al. Effectiveness of XP-endo Finisher and XP-endo Finisher R in removing root filling remnants: a micro-CT study. *Int Endod J.* 2018;51(1):86–91.
23. Zand V, Mokhtari H, Reyhani M-F, Nahavandizadeh N, Azimi S. Smear layer removal evaluation of different protocol of Bio Race file and XP-endo Finisher file in corporation with EDTA 17% and NaOCl. *J Clin Exp Dent.* 2017 Nov;9(11):e1310–4.
24. Turkyaydin D, Demir E, Basturk FB, Övecoglu HS. Efficacy of XP-Endo Finisher in the Removal of Triple Antibiotic Paste from Immature Root Canals. *J Endod.* 2017 Sep 1;43(9):1528–31.
25. Uygun AD, Gündoğdu EC, Arslan H, Ersoy İ. Efficacy of XP-endo finisher and TRUShape 3D conforming file compared to conventional and ultrasonic irrigation in removing calcium hydroxide. *Aust Endod J J Aust Soc Endodontology Inc.* 2017 Aug;43(2):89–93.
26. Wigler R, Dvir R, Weisman A, Matalon S, Kfir A. Efficacy of XP-endo finisher files in the removal of calcium hydroxide paste from artificial standardized grooves in the apical third of oval root canals. *Int Endod J.* 2017;50(7):700–5.
27. Hamdan R, Michetti J, Pinchon D, Diemer F, Georgelin-Gurgel M. The XP-Endo Finisher for the removal of calcium hydroxide paste from root canals and from the apical third. *J Clin Exp Dent.* 2017 Jul;9(7):e855–60.
28. Gokturk H, Ozkocak I, Buyukgebiz F, Demir O. Effectiveness of various irrigation protocols for the removal of calcium hydroxide from artificial standardized grooves. *J Appl Oral Sci Rev FOB.* 2017 Jun;25(3):290–8.
29. Leoni GB, Versiani MA, Silva-Sousa YT, Bruniera JFB, Pécora JD, Sousa-Neto MD. Ex vivo evaluation of four final irrigation protocols on the removal of hard-tissue debris from the mesial root canal system of mandibular first molars. *Int Endod J.* 2017;50(4):398–406.
30. Karamifar K, Mehrasa N, Pardis P, Saghiri MA. Cleanliness of Canal Walls following Gutta-Percha Removal with Hand Files, RaCe and RaCe plus XP-Endo Finisher Instruments: A Photographic in Vitro Analysis. *Iran Endod J.* 2017;12(2):242–7.
31. Elnaghy AM, Mandorah A, Elsaka SE. Effectiveness of XP-endo Finisher, EndoActivator, and File agitation on debris and smear layer removal in curved root canals: a comparative study. *Odontology.* 2017 Apr;105(2):178–83.
32. Bao P, Shen Y, Lin J, Haapasalo M. In Vitro Efficacy of XP-endo Finisher with 2 Different Protocols on Biofilm Removal from Apical Root Canals. *J Endod.* 2017 Feb 1;43(2):321–5.
33. Keskin C, Sarıyılmaz E, Sarıyılmaz Ö. Efficacy of XP-endo Finisher File in Removing Calcium Hydroxide from Simulated Internal Resorption Cavity. *J Endod.* 2017 Jan 1;43(1):126–30.

34. Gokturk H, Ozkocak I, Buyukgebiz F, Demir O. An in vitro evaluation of various irrigation techniques for the removal of double antibiotic paste from root canal surfaces. *J Appl Oral Sci Rev FOB*. 2016 Dec;24(6):568–74.
35. Alves FRF, Andrade-Junior CV, Marceliano-Alves MF, Pérez AR, Rôças IN, Versiani MA, et al. Adjunctive Steps for Disinfection of the Mandibular Molar Root Canal System: A Correlative Bacteriologic, Micro-Computed Tomography, and Cryopulverization Approach. *J Endod*. 2016 Nov 1;42(11):1667–72.
36. Alves FRF, Marceliano-Alves MF, Sousa JCN, Silveira SB, Provenzano JC, Siqueira JF. Removal of Root Canal Fillings in Curved Canals Using Either Reciprocating Single- or Rotary Multi-instrument Systems and a Supplementary Step with the XP-Endo Finisher. *J Endod*. 2016 Jul 1;42(7):1114–9.
37. Azim AA, Aksel H, Zhuang T, Mashtare T, Babu JP, Huang GT-J. Efficacy of 4 Irrigation Protocols in Killing Bacteria Colonized in Dentinal Tubules Examined by a Novel Confocal Laser Scanning Microscope Analysis. *J Endod*. 2016 Jun 1;42(6):928–34.

