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Golden rules

Speed : 600 rpm
Torque : 1.5 Ncm

- Long back and forth gentle strokes.
- Light touch, let the instrument work.
- Work 3-4 seconds in a row, get out, clean the blade and irrigate the canal.

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iRace 1 blister includes / zawiera / περιλαμβάνει: 1 × R1, 1 × R2, 1 × R3

For most cases

- straight, slightly curved and/or large
Prepare glide path first, then

1. Introduce R1 in rotation and reach working length (WL); **In case R1 does not reach the WL, do not force and go to step 1 of iRace Plus protocol**
2. Continue shaping with R2 up to WL;
3. Finish the shaping with R3 up to WL

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iRace Plus 1 blister includes / zawiera / περιλαμβάνει: 2 × (R1a + R1b)

For difficult cases

- severely curved, narrow and/or calcified
Following of iRace protocol step 1

1. Use R1a to reach WL;
 2. Continue shaping with R1b up to WL;
- When WL is reached, continue from step 2 of iRace protocol.

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EN Sterilization protocol for reusable medical devices

Non sterile products

First use

Carefully follow the rules herebelow
Cleaning and Sterilization

Sterile products

First use

1. Open the alveolus
2. Extract the instrument (use of gloves imperative)
3. Use the instrument according to good dental practice

Second use:

Carefully follow the rules herebelow
Second and following uses.

Cleaning

1. Manual cleaning with/without Ultrasonic assistance

Thorough cleaning allows efficient disinfection & sterilization.

2. Rinse

Under distilled/demineralized water at least 1 min.
Then dry the instruments.

3. Inspection - Check
Sort out damaged or worn out instruments.

4. Packing

Place the instruments in a suitable support/ container and pack the devices in sterilisation pouches ISO 11607-1.

Sterilization

5. According EN/ISO 17664 protocol

Autoclave: 134°C/273°F,
2,2 bar during 18 min.

6. Storage

Keep devices in sterilization wrap/pouch in a dry and clean environment.

 Always refer to manufacturer's instructions for cycle & duration

Second and following uses

7. Pre-disinfection

Plunge instruments after use in detergent solution and eventually brush them manually.

8. Rinse

Under running water at least 1 min.
Then dry the instruments.

9. Follow steps 1 - 6 above

 Complete operating protocol and warnings on www.fkg.ch

General information

Pre-/Disinfection

DO NOT USE

- solutions containing : phenol (corrosion)
- aldehyde (blood fixation)
- di-/trietanolamines (corrosion)

Sterilization with chemiclav or dry air devices have not been validated by FKG Dentaire.
Refer to manufacturer's instructions for cycle & duration. Nevertheless our instruments withstand such methods.

Caution:

This product contains nickel and should not be used for individuals with known allergic sensitivity to this metal.

A single patient use is recommended to avoid cross-contamination.

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