

Speed & torque

Product	Speed - minimum	Speed - optimal	Torque
BioRace	600	600	1.5
BioRace Extended	600	600	1.5
BT-Apisafe	600	800	1
BT-Race	600	800	1.5
BT-Race XL	600	800	1.5
D-Race DR1	800	1'000	1.5
D-Race DR2	400	600	1.5
EasyRace	400	600	1.5
iRace	400	600	1.5
iRace Plus	400	600	1.5
PreRace	400	600	1.5
Race	400	600	1.5
Race ISO10	600	800	1.5
ScoutRace	600	800	1.5
XP-endo Finisher	800	1000	1
XP-endo Finisher R	800	1000	1
XP-endo Shaper	800	1000	1
XtremRace	400	600	1.5