



FKG  
swiss endo




### ► BT-Apisafe

To prepare in complete safety the apical third inside severely curved, narrow or sclerotic canals, while preserving the coronal third.




Zero taper  
for greater flexibility



BT-Tip  
for unparalleled efficiency

Availability : ISO 25/30/40/50/60



1639

FKG Dentaire Sàrl  
Le Crêt-du-Locle 4  
1639 Locle  
Switzerland

FKG REF: 996AA.00.02A-EN\_REV.4\_2022-01 - © FKG Dentaire Sàrl - Pictures of our products are not contractual.

FKG Dentaire Sàrl  
www.fkg.ch