

Press release

TotalFill® BC Sealer HiFlow™

FKG Dentaire presents its new product for warm obturation

La Chaux-de-Fonds, Switzerland, March 2019.

The "State-of-the-Art" in endodontic obturation is now available for warm techniques !

Clinicians now have two distinct formulations of the same proven pre-mixed bioceramic sealer to choose from. Compared to standard **TotalFill® BC Sealer™**, **TotalFill® BC Sealer HiFlow™** exhibits a lower viscosity when heated and is more radiopaque, making it optimized for warm obturation techniques.

Top Features:

- Condensation-based technique
- Optimal for warm obturation
- Lower viscosity when heated
- Higher radiopacity



Standard gutta-percha can be used with **TotalFill® BC Sealer™**, but for a true, gap-free seal, **TotalFill® BC Points™** are recommended. **TotalFill® BC Points™** are impregnated and coated with bioceramic particles to allow for bonding with **TotalFill® BC Sealers**. For warm techniques it is



recommended to use **TotalFill® BC Sealer HiFlow™** with the **TotalFill® BC Points™ 150 Series**, along with **TotalFill® BC Pellets™** at 150 °C – 220 °C.

For further information, contact the FKG Dentaire team:

FKG Dentaire SA

Crêt-du-Loche 4

CH-2304 La Chaux-de-Fonds

Switzerland

T. +41 32 924 22 44

F. +41 32 924 22 55

info@fkg.ch

www.fkg.ch

www.facebook.com/FKGDentaire

About FKG Dentaire SA:

Founded in Switzerland in 1931, FKG Dentaire SA gained a new momentum in 1994, when the Rouiller family took over the reins of the company. The family propelled FKG to the forefront in the development, manufacturing and distribution of dental products destined for general practitioners, endodontists and laboratories. The FKG strategy is centered on innovative high-precision products and the creation of machines designed specifically for the dental field. Its aim is to offer solutions that meet the most demanding needs of end users.