



Press release

XP-endo® - 3D Solutions for Root Canal Treatment

FKG Dentaire presents its new 3D Solutions

La Chaux-de-Fonds, Switzerland, March 2019.

Since launching its first revolutionary 3D instrument – the XP-endo® Finisher – four years ago, Swiss manufacturer of innovative, high quality endodontic solutions, FKG Dentaire SA, introduces the latest application expansions to its total therapy system : XP-endo® Solutions.

XP-endo® Solutions comprises 2 different files : XP-endo® Shaper and XP-endo® Finisher, that expand at body temperature to achieve safe and easy management of root canal therapy while preserving dentine.

Together with the introduction at IDS 2019 of XP-endo® Retreatment, a new clinical sequence focused on more efficient retreatment solutions, XP-endo® Solutions delivers a more complete system to both practitioners and patients in initial treatment and retreatment management, requiring fewer files.

A Complete System

XP-endo® Shaper, with its snake-like shape, is primarily dedicated to canal shaping. Thanks to its adaptive core and a 1% taper associated to a guiding tip, XP-endo® Shaper is able to shape all types of canals with greater safeguarding of the original shape of the root – from regular structures to the most complex. The XP-endo® Retreatment clinical sequence utilizes XP-endo® Shaper that acts like a corkscrew to remove filling material from the canal with extreme precision and efficiency.

XP-endo® Finisher is designed to provide optimal expansion properties through a minimalistic design that ensures smooth and extraordinary, efficient cleaning capabilities of root canal walls.



XP-endo® Finisher is available in 2 sizes : ISO 025 and ISO 030 (Finisher R), both with a 0% taper. The 2 instruments can clean canal walls, with the Finisher R (Reinforced) dedicated to retreatment cases and the smaller version (Finisher) intended for initial treatment use.

XP-endo® Solutions is by far the most complete system available for 3D endodontics today offering highly adaptive, less invasive treatment solutions with less files... *Endo Done !*

www.fkg.ch/xpendo

For further information, contact the FKG Dentaire team:

FKG Dentaire SA

Crêt-du-Loche 4

CH-2304 La Chaux-de-Fonds

Switzerland

T. +41 32 924 22 44

F. +41 32 924 22 55

info@fkg.ch

www.fkg.ch

www.facebook.com/FKGDentaire

About FKG Dentaire SA:

Founded in Switzerland in 1931, FKG Dentaire SA gained a new momentum in 1994, when the Rouiller family took over the reins of the company. The family propelled FKG to the forefront in the development, manufacturing and distribution of dental products destined for general practitioners, endodontists and laboratories. The FKG strategy is centered on innovative high-precision products and the creation of machines designed specifically for the dental field. Its aim is to offer solutions that meet the most demanding needs of end users.