


**FKG Dentaire introduces 3 new rotary NiTi instruments:
Glide path files RaCe ISO 10 .02/.04/.06**

Description

Thanks to their extremely small tip RaCe ISO 10 can be used directly after manual instruments ISO 08. Considerable time saving in all kind of canals!

The RaCe ISO 10 offer the following features:

- Exclusive rounded Safety Tip
- Unique anti-screwing design
- Sharp cutting edges for efficiency
- Electro-chemical polishing: better resistance to torsion & fatigue
- *SafetyMemoDisc*: for mastering metal fatigue & controlling number of uses



RaCe : the most complete range of rotary NiTi endodontic instruments! From mechanical glide path to final shaping...

		Ø ISO										
		Taper SMD										
		10	15	20	25	30	35	40	45	50	55	60
PRE-RaCe 19 mm	.10							●				
	.08						●					
	.06					●		●				
RaCe 25 mm (21/31 mm)	.06	●	●	●	●	●	●	●				
	.04	●	●	●	●	●	●	●		●		
	.02	●	●	●	●	●	●	●	●	●	●	●

Availability

- Individual sizes (5 pces/box)